---

# ZIMMERMANN

---

## IMMOB-SCHIENE

Immob Splint

comfort

## Sprachen

Languages

<b>DE</b> deutsch	Gebrauchsanweisung .....	02
<b>EN</b> english	Instructions for use .....	05

## **DE** Immob-Schiene

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## **EN** Immob Splint

Many thanks for placing your trust in this medical device. Please read the existing instructions for use carefully. If you have any questions, please refer to your physician or the specialist retailer from whom you purchased this medical device.

### **Intended purpose**

This medical device is a brace for immobilisation of the knee joint in a defined position (bent/stretched).

### **Indications**

Pre-operative, post-operative, post-traumatic, joint, muscle or ligament injuries

### **Contraindications**

Peripheral arterial occlusive disease (PAOD), lymph drainage disorders, also unclear soft tissue swellings distal to the aid positioned, sensory and circulatory disorders with the region of the body treated, skin diseases in the part of the body treated.

### **Application risks/Important notes**

This medical device is a prescribable product. Discuss the use and duration with your treating physician. The expert staff from whom you have received the medical device will select the appropriate size and instruct you regarding its use.

- remove the medical device prior to radiological examinations
- in case of numbness or pressure sores, loosen the medical device and remove it if necessary
- in case of persistent complaints, consult the physician or specialist retailer
- use the medical device in accordance with therapeutic needs
- only use other products simultaneously after consultation with your physician

- do not make any changes to the product
- do not wear it on open wounds
- do not use in case of intolerance of one of the materials used
- no re-use – this medical aid is intended for treating one patient

### **Fitting instructions**

The Immob Splint with Patella Recess consists of 3 parts: A large middle part and 2 side parts attached with Velcro.

Before attaching the brace, the integrated aluminium splints which are in the splint bags for the middle and side parts, must be formed to fit the patient's existing body contours. The desired bending knee joint position can be re-adjusted by forming the aluminium splints on the middle part of the brace. The freely positionable side parts enable adaptation to various circumference measurements. When set correctly, the side parts should overlap by approx. 2 cm (0.79 inches). Lay the Immob Splint around the leg from posterior so that the kneecap (patella) is enclosed by the hollow. Thread Velcro fasteners through plastic eyelets, guide them back and attach them with Velcro under slight tension. Attach the Immob Splint so that, although the knee joint is immobilised, no congestions are caused in the calf and foot area.

### **Removal**

To remove the product, loosen the belt straps and take the brace off the leg.

### **Material composition**

Polyamide (PA), cotton (CO), PUR foam (PUR)

The sewn in textile label provides the precise material composition.

### **Cleaning information**

 Delicate wash  Do not clean chemically  Do not bleach

 Do not dry in a tumble dryer  Do not iron

Do not use fabric conditioner. Close the Velcro fastener in order to avoid damaging other items of clothing and remove the aluminium splints. Stretch back into shape and dry in the fresh air.

### **Guarantee**

The legal regulations of the country in which you acquired the product apply to the purchased product. Please contact your

specialist retailer if you suspect a warranty claim. Please clean the product before submitting a warranty claim. If the enclosed instructions for use have not been properly observed, the warranty may be impaired or cancelled. The warranty does not cover use of the product inappropriate for the indication, non-observance of application risks, instructions and unauthorised modifications to the product.

### **Useful life/Lifetime of the product**

The lifetime of the medical device is determined by its natural wear and tear if treated appropriately and as recommended.

### **Duty of notification**

If a serious worsening of the state of health occurs when using the medical device, you are obliged to notify the specialist dealer or the manufacturer and the Medicines & Healthcare products Regulatory Agency (MHRA). You can find the manufacturer's contact data in these instructions for use. You can find the contact information for the appointed authority for your country under the following link: [www.bort.com/md-eu-contact](http://www.bort.com/md-eu-contact).

### **Disposal**

Upon the termination of use, the product must be disposed of in accordance with the corresponding local requirements.

### **Declaration of conformity**

We confirm that this device conforms with the requirements of REGULATION (EC) NO. 2017/745 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL.

Status: 10.2019

Medical device |  Single patient – multiple use

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**Zimmermann Sanitäts- und Orthopädiehaus GmbH**
**BORT GmbH**